NGS Data Analysis

Roberto Preste

Useful info

Contacts:

roberto.preste@gmail.com


Slides:

http://bit.ly/ngs-data

NGS data analysis

Overview

NGS Data Analysis: the basic idea




NGS data

R1

R3

Organization and analysis of NGS data

NGS Data Analysis: the actual workflow



NGS data analysis

Quality check & preprocessing

Fasta files

>J01415.2 Homo sapiens mitochondrion, complete genome GATCACAGGTCTATCACCCTATTAACCACTCACGGGAGCTCTCCATGCATTTGGTATTTTCGTCTGGGGG GTATGCACGCGATAGCATTGCGAGACGCTGGAGCCCGGAGCACCCTATGTCGCAGTATCTGTCTTTGATTC ATTAATTAATGCTTGTAGGACATAATAATAACAATTGAATGTCTGCACAGCCACTTTCCACACAGACATC ATAACAAAAATTTCCACCAAACCCCCCCCCCCCCGCTTCTGGCCACAGCACTTAAACACATCTCTGCCA

Sequence ID

Sequence

- Both human- and machine-readable
- Can store multiple sequences
- ID can contain details or comments
- Usually contains full genomes or long sequence chunks

Fastq files

- Both human- and machine-readable
- Can store multiple sequences
- ID can contain sequencing details and technical info
- Usually contains short sequence chunks (sequencing reads)

Quality score

Phred quality score: estimated probability of an error in base calling

usually [0-40] —
$$Q_{sanger} = -10 \log_{10} P$$
 — Probability that the base call is incorrect

Encoded using ASCII characters in fastq files:

Quality score	Probability of errors	ASCII encoding
0-9	1	!"#\$%&'()*
10-19	1/10	+,/01234
20-29	1/100	56789:;<=>
30-39	1/1000	?@ABCDEFGH
40	1/10000	I

Quality score

@HWUSI-EAS100R:6:73:941:1973#0/1

GATTTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT

+

!''*((((***+))%%%++)(%%%).1***-+*''))**55CCF>>>>>CCCCC420

0 6 9

Quality score	Probability of errors	ASCII encoding
0-9	1	!"#\$%&'()*
10-19	1/10	+,/01234
20-29	1/100	56789:;<=>
30-39	1/1000	?@ABCDEFGH
40	1/10000	l l



Quality check

FastQC: visual report of several quality checks for NGS data, useful for further processing




FastQC modules


Basic Statistics

Measure	Value
Filename	ionXpress002.fastq
File type	Conventional base calls
Encoding	Sanger / Illumina 1.9
Total Sequences	65578765
Sequences flagged as poor quality	0
Sequence length	25-392
%GC	52

Sequence Length Distribution







Cleansing of reads to solve several issues:

remove adapters



Cleansing of reads to solve several issues:


Good quality

Bad quality

Adapter

remove adapters

cut low quality bases from both ends



Cleansing of reads to solve several issues:

Good quality
Bad quality
Adapter

remove adapters

cut low quality bases from both ends

drop short reads



Cleansing of reads to solve several issues:

remove adapters

cut low quality bases from both ends

drop short reads

Common tools:

- Trimmomatic
- Trim Galore!
- FASTX


Post-processing quality check



Measure	Value
Filename	ionXpress002.trimmed.fast
File type	Conventional base calls
Encoding	Sanger / Illumina 1.9
Total Sequences	34710199
Sequences flagged as poor quality	0
Sequence length	30-307
%GC	53

Was the processing effective?

Are these data ready to be aligned?



NGS data analysis

Alignment

Alignment vs Assembly

Alignment (reference-based)

reference genome available

→ reads aligned on it

Assembly (de-novo)

reference genome not available

GTGACTTAGT

→ reads aligned with each other

GCTAGCTAGT

..GTGACTTAGTCGTAGCTAGCTAGTAGCTCGATCTAGA..

GTGACTTAGT TAGTAGCTCG

GAGTTAGTCG GTAGCTCGAT

CTTAGTCGTA AGCTCGACCT

TAGTCGTAGC CTCGACCTAG

AGTTAGTCGT GTAGCTCGAT

TCGTAGCTAG AGCTCGACCT

TGAGTTAGCC AGCTAGTAGC

..GTGACTTAGTCGTAGCTAGCTAGTAGCTCGATCTAGA..

CGATCTAGA

Alignment

..GTGACTTAGTCGTAGCTAGCTAGTAGCTCGATCTAGA..

Reference genome

GTGACTTAGT

TAGTAGCTCG

GAGTTAGTCG

GTAGCTCGAT

CTTAGTCGTA

AGCTCGACCT

TAGTCGTAGC

CTCGACCTAG

Common aligners:

- BWA
- Bowtie
- GMAP/GSNAP





http://bit.ly/2DG8IQj

Assembly

http://bit.ly/2DnwRdt

GTGACTTAGT GCTAGCTAGT CGATCTAGA

AGTTAGTCGT GTAGCTCGAT

TCGTAGCTAG AGCTCGACCT

TGAGTTAGCC AGCTAGTAGC

Reads

..GTGACTTAGTCGTAGCTAGCTAGTAGCTCGATCTAGA..

Common approaches:

- greedy algorithm
- graph method

Consensus sequence

Common assemblers:


- Newbler
- SPAdes
- MaSuRCA

Paired-end vs single-end reads

Easy to identify read position in genome

High accuracy for structural rearrangements and assembly of repetitive regions

More expensive



Suitable for most applications

Cheaper

Faster

SAM/BAM files

Text-based format used to store aligned reads (to a reference genome)

```
QHD VN:1.6 SO:coordinate
                                    39 TTAGATAAAGGATACTG *
r002
                                 O O AAAAGATAAGGATA
       0 ref 9 30 3S6M1P1I4M *
r003
      0 ref 9 30 5S6M
                                 O O GCCTAAGCTAA
                                                         * SA:Z:ref,29,-,6H5M,17,0;
r004
       0 ref 16 30 6M14N5M
                                     O ATAGCTTCAGC
r003 2064 ref 29 17 6H5M
                                                         * SA:Z:ref,9,+,5S6M,30,1;
                                     O TAGGC
     147 ref 37 30 9M
                                   -39 CAGCGGCAT
                                                         * NM:i:1
```


SAM (Sequence Alignment Map)

- Both human- and machine-readable
- Header section contains TAG:VALUE pairs
- Alignment section contains 11 mandatory fields

BAM (Binary Alignment Map)

- Compressed version of SAM
- Binary format
- Only machine-readable

SAM files



@HD: header line delimiter

VN: version number

S0: alignments sorting order

@SQ: reference sequence

SN: reference sequence name

LN: reference sequence length


SAM files

QNAME	Ē		N	MAP	Q T	SEQ					
r001	99	ref	7	30	8M2I4M1D3M	=	37	39	TTAGATAAAGGATACTG	*	
r002	0	ref	9	30	3S6M1P1I4M	*	0	0	AAAAGATAAGGATA	*	
r003	0	ref	9	30	5S6M	*	0	0	GCCTAAGCTAA	*	SA:Z:ref,29,-,6H5M,17,0;
r004	0	ref	16	30	6M14N5M	*	0	0	ATAGCTTCAGC	*	
r003	2064	ref	29	17	6H5M	*	0	0	TAGGC	*	SA:Z:ref,9,+,5S6M,30,1;
r001	147	ref	37	30	9M	=	7	-39	CAGCGGCAT	*	NM:i:1

POS	CIGAR

Col	Field	Type	Regexp/Range	Brief description
1	QNAME	String	[!-?A-~]{1,254}	Query template NAME
2	FLAG	\mathbf{Int}	$[0, 2^{16} - 1]$	bitwise FLAG
3	RNAME	String	* [!-()+-<>-~][!-~]*	Reference sequence NAME
4	POS	\mathbf{Int}	$[0, 2^{31} - 1]$	1-based leftmost mapping POSition
5	MAPQ	Int	$[0, 2^8 - 1]$	MAPping Quality
6	CIGAR	String	* ([0-9]+[MIDNSHPX=])+	CIGAR string
7	RNEXT	String	* = [!-()+-<>-~][!-~]*	Ref. name of the mate/next read
8	PNEXT	Int	$[0, 2^{31} - 1]$	Position of the mate/next read
9	TLEN	Int	$[-2^{31}+1, 2^{31}-1]$	observed Template LENgth
10	SEQ	String	* [A-Za-z=.]+	segment SEQuence
11	QUAL	String	[!-~]+	ASCII of Phred-scaled base QUALity+33

Alignment quality check



Integrative Genomics Viewer (IGV)

Interactive visualization of NGS data from Fasta, SAM, BAM, VCF files


Additional features are organized in tracks (gene expression, methylation, copy number variations...)

NGS data analysis

Variant calling


Variations

http://bit.ly/2QHQocd



Aligned data can be used to assess the presence of mutations:

- somatic
- germline



Variations

Most common variations searched for:

\/_	⊬: ~ +	:	4	
Va	Πaι	ion	ιv	рe
			- 1	


Reference ACTGACGCATGCATGCATGC

SNP ACTGACGCATGCATCATTCATGC

Insertion ACTGACGCATGGTACATGCATGC

Deletion ACTGACGC--GCATCATGCATGC

Variation effect



http://bit.ly/2K9a3zj

Variant callers

A plethora of different tools, each with its own peculiarities:

- variation type (SNV vs structural variation)
- source (somatic vs cancer mutations)
- only detect variants vs also predict their effect

















Variant Call Format (VCF)

VCF files

http://bit.ly/2Kfl80W

```
##fileformat=VCFv4.2
##fileDate=20090805
##source=myImputationProgramV3.1
##reference=file:///seq/references/1000GenomesPilot-NCBI36.fasta
##contig=<ID=20,length=62435964,assembly=B36,md5=f126cdf8a6e0c7f379d618ff66beb2da,species="Homo sapiens",taxonomy=x>
##phasing=partial
##INFO=<ID=NS, Number=1, Type=Integer, Description="Number of Samples With Data">
##INFO=<ID=DP, Number=1, Type=Integer, Description="Total Depth">
##INFO=<ID=AF, Number=A, Type=Float, Description="Allele Frequency">
##INFO=<ID=AA, Number=1, Type=String, Description="Ancestral Allele">
##INFO=<ID=DB, Number=0, Type=Flag, Description="dbSNP membership, build 129">
##INFO=<ID=H2,Number=0,Type=Flag,Description="HapMap2 membership">
                                                                                       Header
##FILTER=<ID=q10,Description="Quality below 10">
##FILTER=<ID=s50,Description="Less than 50% of samples have data">
                                                                         (information and metadata)
##FORMAT=<ID=GT, Number=1, Type=String, Description="Genotype">
##FORMAT=<ID=GQ, Number=1, Type=Integer, Description="Genotype Quality">
##FORMAT=<ID=DP, Number=1, Type=Integer, Description="Read Depth">
##FORMAT=<ID=HQ.Number=2.Type=Integer.Description="Haplotype Quality">
#CHROM POS
               TD
                         REF
                                ALT
                                        QUAL FILTER INFO
                                                                                      FORMAT
                                                                                                   NA00001
                                                                                                                  NA00002
20
       14370
              rs6054257 G
                                        29
                                             PASS
                                                    NS=3:DP=14:AF=0.5:DB:H2
                                                                                      GT:GQ:DP:HQ 0|0:48:1:51,51 1|0:48:8:51,51
       17330
                                             q10
20
                                                    NS=3;DP=11;AF=0.017
                                                                                      GT:GQ:DP:HQ 0|0:49:3:58,50 0|1:3:5:65,3
                                G,T
                                                    NS=2;DP=10;AF=0.333,0.667;AA=T;DB GT:GQ:DP:HQ 1|2:21:6:23,27 2|1:2:0:18,2
20
       1110696 rs6040355 A
                                             PASS
                                             PASS
                                                    NS=3:DP=13:AA=T
                                                                                      GT:GQ:DP:HQ 0|0:54:7:56.60 0|0:48:4:51.51
20
       1230237 .
20
       1234567 microsat1 GTC
                                G, GTCT
                                        50
                                             PASS
                                                    NS=3;DP=9;AA=G
                                                                                      GT:GQ:DP
                                                                                                  0/1:35:4
                                                                                                                  0/2:17:2
```

Variants

Variant annotations

Genotype annotations

VCF files

Mandatory fields:

#CHROM: chromosome or contig

POS: variant position (1-based)

ID: variant identifier (usually dbSNP ID)

REF: reference base(s)

Optional fields:

FORMAT: genotype information fields

ALT: alternative base(s)

QUAL: Phred quality score for each ALT

FILTER: variant call filter status

INFO: key-value pairs with additional information (described in header)

SAMPLE1 ... SAMPLEn: values for fields listed in FORMAT

NGS data analysis

Functional annotation

Functional annotation

~3.000.000 small variants (SNP/indels) per genome

ATCATGCATGC

ATCATTCATGC

Which ones are most interesting for our purpose?







http://bit.ly/2QekbN3



http://bit.ly/2BfefuO



Annovar

Identify variants and flag those with detrimental effects

ANNOVAR: functional annotation of genetic variants from high-throughput sequencing data

Kai Wang¹,*, Mingyao Li² and Hakon Hakonarson^{1,3}

http://bit.ly/2zc9fp5

Gene-based annotations

- → identify variants that can cause protein coding changes
- → detect the amino acids that are affected

Annovar

Identify variants and flag those with detrimental effects

ANNOVAR: functional annotation of genetic variants from high-throughput sequencing data

Kai Wang^{1,*}, Mingyao Li² and Hakon Hakonarson^{1,3}

http://bit.ly/2zc9fp5

Region-based annotations

- → identify variants in specific genomic regions:
 - conserved regions
 - predicted transcription factor binding sites
- segmental duplication regions
- etc

Annovar

Identify variants and flag those with detrimental effects

ANNOVAR: functional annotation of genetic variants from high-throughput sequencing data

Kai Wang¹,*, Mingyao Li² and Hakon Hakonarson^{1,3}

http://bit.ly/2zc9fp5

Filter-based annotations

- → identify variants that are documented in specific databases:
 - dbSNP
 - 1000 Genome Project
 - etc

Annovar VCF files

Extensive number of new annotations added to the initial VCF file

Chr	Start	t	End	Ref	Alt	Func.refGe	Gene.refGene	GeneDetail.refGene	ExonicFunc.refGene	AAChange.refGene	Xref.refGene	ExAC_Freq
	1 9	948921	948921	Т	С	UTR5	ISG15	NM_005101:c33T>C			Immunodeficienc	0.941
	1 14	404001	1404001	G	T	UTR3	ATAD3C	NM_001039211:c.*91G>T				0.054
	1 59	935162	5935162	Α	T	splicing	NPHP4	NM_001291594:exon17:c.1282-2T>A		20	Nephronophthisi	0.825
	1 1627	736463	162736463	C	T	intronic	DDR2	AND THE PROPERTY OF THE PARTY O			Spondylometaepi	
	1 848	875173	84875173	С	T	intronic	DNASE2B	2				2
	1 132	211293	13211294	TC	-	intergenic	PRAMEF36P;F	dist=11566;dist=116902				
	1 114	403596	11403596	-	AT	intergenic	UBIAD1;PTCH	dist=55105;dist=135699				
	1 1054	492231	105492231	Α	ATAAA	intergenic	LOC10012913	dist=872538;dist=640085				
	1 677	705958	67705958	G	A	exonic	IL23R		nonsynonymous SNV	IL23R:NM_144701:exon9:c.G1142A:p.R381Q	*	0.041
	2 2341	183368	234183368	Α	G	exonic	ATG16L1		nonsynonymous SNV	ATG16L1:NM_198890:exon5:c.A409G:p.T137A;ATG16L1:NM_0:		0.457
1	6 507	745926	50745926	С	T	exonic	NOD2		nonsynonymous SNV	NOD2:NM_001293557:exon3:c.C2023T:p.R675W;NOD2:NM_0	Blau syndrome, A	0.023
1	6 507	756540	50756540	G	C	exonic	NOD2	2	nonsynonymous SNV	NOD2:NM_001293557:exon7:c.G2641C:p.G881R;NOD2:NM_0:	Blau syndrome, A	0.009917
1	6 507	763778	50763778	-3	С	exonic	NOD2	2	frameshift insertion	NOD2:NM_001293557:exon10:c.2936dupC:p.L980Pfs*2;NOD2	Blau syndrome, A	0.013
1	3 207	763686	20763686	G	-	exonic	GJB2		frameshift deletion	GJB2:NM_004004:exon2:c.35delG:p.G12Vfs*2	Bart-Pumphrey s	0.006038
1	3 207	797176	21105944	0	-	exonic	CRYL1;GJB6		frameshift deletion	GJB6:NM_001110220:wholegene;GJB6:NM_001110221:whole		
	8 88	887543	8887543	Α	Т	exonic	ERI1		stoploss	ERI1:NM_153332:exon7:c.A1049T:p.X350L	*	,
	8 88	887539	8887539	Α	T	exonic	ERI1	*	stopgain	ERI1:NM_153332:exon7:c.A1045T:p.K349X	89	
	8 88	887536	8887537	AG	GATT	exonic	ERI1		frameshift substitution	ERI1:NM_153332:exon7:c.1042_1043GATT:p.R348Dfs*2	37	9
	8 88	887540	8887540	G	GGAA	exonic	ERI1	÷	nonframeshift substit	ERI1:NM_153332:exon7:c.1046delinsGGAA:p.R348_K349insR		
	5 12	295288	1295288	G	A	upstream	TERT	dist=126			31	gr


http://bit.ly/2Q6cvfr

MToolBox

Human mtDNA reconstruction, analysis and annotation from NGS data

Haplogroup predictions

Both command-line and web-based versions available





HmtDB

Over 48000 human mitochondrial sequences

Healthy/patient and continent-specific subsets

Genome-centric





Human Mitochondrial DataBase


https://www.hmtdb.uniba.it

HmtVar

Over 40000 human mitochondrial variants

Pathogenicity predictions available for most variants

Variant-centric



Main information regarding the T3440G variant.

Variant Card

Main Info

Main Info

Basic

Variability and pathogenicity data available for T3440G.

Pathogenicity Predictions External Resources

NGS Data Analysis: pipelines



General-purpose programming language

Easy to learn

Very powerful (libraries available for anything you can think of)

→ Biopython: specific module for bioinformatics

NGS Data Analysis: pipelines

Originally suited for statistics

Particularly used for data analysis and visualization

Gained a lot of traction for many different disciplines

→ Bioconductor: specific package for bioinformatics



Useful info

Contacts:

roberto.preste@gmail.com

Slides:

http://bit.ly/ngs-data